

**HPLC SOP**

Door code: 1-34 (as of December, 2008)  
Username: 3000Hanover  
Password: 3000Hanover

1. Turn on all Hardware (instrument buttons, 6 in total) before the software (computer)
  - a. Start at top and work down and right
  - b. Hear clicking noises that indicates the start of the pump cycle
  - c. To end, make sure Hardware is OFF before the Software is turned OFF
2. Click on the icon for Software Instrument 1 on-line
  - a. If software does not load or there is an error, log out and restart the computer
3. Software shows syringe-pump-column-detector, all can be clicked on to change and monitored
4. Sample icon is top left of screen, single or multiple see 1 bottle or numerous bottles
5. Mobile Phase
  - a. Individual containers for each solvent
  - b. Follow the lines (numbered) from the bottles to the inlet into the column
  - c. Milli-Q water only!
  - d. MeOH
  - e. Make sure containers are full and the software knows how much volume in the containers
  - f. If solvent is less than 200 mL, no more injections will be performed
6. Check for waste flow either
  - a. Empties into the waste container, make sure tubing is in the bottle and not the floor
  - b. Flows into 2<sup>nd</sup> detector for further analysis
7. Moisture Check-there is an automatic shut-off if the instrument gets too wet.
  - a. Always check for leaks and excess moisture at every fitting
8. Flush the column for 20-30 mins using strong to weak solvent
  - a. 100 % H<sub>2</sub>O
  - b. 100 % MeOH
  - c. mixture of H<sub>2</sub>O-MeOH
9. Software Pump, click on icon
  - a. 1 mL/min flow rate
  - b. Increase flow very slowly to not increase the pressure on the column to rapidly, do NOT EXCEED 400 bar maximum pressure
  - c. Stoptime = 30 min
  - d. Control, click on OK, hear pump turn on
10. Detector, click on icon, depending on matrices and background of sample
  - a. Data is stored only to the detector that has the box checked
  - b. All data is stored to the end of the run
  - c. Stop Time 40 min
  - d. Diode Array detector: Range 210-400 nm
    - i. Default setting are OK for running
  - e. ADC is always ON, not saved, data not needed

- i. Set-up peaks for integration (set-up will prevent integrating every wiggle)
  - ii. Width 100
  - iii. Reject 1000
  - iv. XXX 1000
11. Column Thermostat
  - a. 30°C is good for a little heat, will feel heat sink warm up that holds the column
  - b. 50-60°C is maximum temperature, if needed
12. Injection Volume
  - a. Smaller is better for better resolution, separation, and sensitivity
  - b. 3-100 µL is range
  - c. Start with 10 µL then try 5 µL to 3 µL is concentrated
  - d. 25 µL is sample is dilute
13. Always make sure solvent is moving through the lines
  - a. Always check for bubbles or air pockets in the sample lines!
  - b. If bubble/no solvent flow, OPEN the by-pass valve, black knob, 2 turns to remove all pressure from the system
  - c. Increase pump speed to 5 mL/min to flush out the lines
  - d. Makes sure to return the pump speed back to 1 mL/min before starting again!
  - e. Close by-pass valve to re-establish backpressure
14. Samples
  - a. Degasses samples are the best, less bubbles
  - b. Make sure all samples are filtered before sending through the column, the column is a very expensive filter
  - c. Use a syringe and ~25 µm filter
15. Save Method
  - a. Save as System File
  - b. Load method
    - i. File-load-method
  - c. Copy spectrum to clipboard, paste into paint, and save a jpg